

TRADITIONAL “PORK&APPLE” SAUSAGE

GLUTEN FREE

Overview

Popular sausage variety, made with CMT gluten free crumbs. Lots of apples. Easy to make and always consistent.

Based on CMT’s revolutionary gluten free “pinhead rusk”

Materials

1. MEAT		
Pork 90 VL, shoulder		4.25 kg
Pork 50 VL, belly not too fat		<u>3.25 kg</u>
		8.00 kg
2. WATER/ICE		1.50 kg
3. INGREDIENTS		
CMT gluten free Cumberland Sausage Mix 600-SM-1412		1.36 kg
	Batch total	11.00 kg approx
4. CMT hog casings 30/32 or 32/35		

Process

1. Mince meat through 5 mm plate. If a coarser cut is desired, then use 7 mm or 10 mm plate.
2. Add seasoning mix (1 bag for 11 kg mix) and water and mix very well until the water is absorbed. Please note that the crumbs do not take up water as quickly as rusk.
3. Remince the mix through 3 mm or 5 mm plate (optional)
4. Fill into natural casings

Useful Tips

1. Always use well-chilled meat and water and keep mincer knives and plates sharp.
2. Always store sausage refrigerated

The above recipe and processing information is accurate to the best of our knowledge. Yet no responsibility or liability can be accepted for any eventual product failure

CMT Continental Meat Technology, 31 Salford Road, Aspley Guise, Milton Keynes MK17 8HT
 Tel: 01908 584489 Fax: 01908 584317 Email: juergenmaurer@c-m-t.demon.co.uk